


# moisture measurement and control...

**CONTINUOUSLY-  
AUTOMATICALLY**

Model W, Electrolytic  
Moisture Analyzer,  
Type NEP-PL

Lightweight, portable moisture analyzer complete with all necessary components for accurate, precise moisture measurement.



**COMPACT** . . . Small in size, 9" wide, 8" deep, 17" high, weight—24 lbs. Portable, yet ruggedly constructed for long, dependable service. (Other types available in panel mounted or explosion-proof housings.)

**ACCURATE** . . . Assures the most accurate readings possible, for ranges 10, 30, 100, 300 and 1000 ppm H<sub>2</sub>O v/v. Equipped with precision sample flow and by-pass controls and integral sample and by-pass flowmeters.

**VERSATILE** . . . Measures and controls moisture in gas streams including manufactured, natural, refrigerant gas, etc. For complete information, write for Data Sheet DS-1010.

**Manufacturers Engineering  
and Equipment Corp.**

Dept. 53-AC 250 Titus Ave., Warrington, Pa.


Circle No. 120 on Readers' Service Card

**New! BY MISCO**

## MICRO CELL HOLDER AND MICRO CELLS

FOR THE BAUSCH & LOMB SPECTRONIC  
20 PHOTO COLORIMETER


NOW AVAILABLE: Adapt your Spectronic 20 to micro samples with the new MISCO micro cell holder and micro cells. Three types of micro cells are available. Determinations can be made with as little as 200 micro-liters total volume. The micro cell can be emptied automatically at the touch of a button, allowing measurement of samples with speed and precision without removing the cell itself.



**REPLACE IT AS A UNIT!**

Get more use from your  
Spectronic 20

Go to the MISCO Micro  
System




**MICROCHEMICAL SPECIALTIES CO.**  
1825 EASTSHORE HIGHWAY BERKELEY 10, CALIFORNIA

Circle No. 185 on Readers' Service Card

QUARTZ PRISM

## MONOCHROMATOR FOR THE VISIBLE AND ULTRAVIOLET



Model  
QPM 30

Miniature monochromator with uni-directional positioned entrance and exit optics for direct insertion into existing light beams and systems.

WAVELENGTH RANGE  
200-700 millimicrons  
SLITS  
Adjustable in unison  
EFFECTIVE APERTURE  
F/7.6

**SCHOEFFEL  
INSTRUMENT CO.**  
Hillsdale, New Jersey

MANUFACTURER OF  
INSTRUMENTS AND  
COMPONENTS FOR THE  
PHYSICAL SCIENCES

Circle No. 153 on Readers' Service Card

## Applications Specialist in Mass Spectrometry

Excellent position for individual with a minimum of a BS degree, advanced degree preferred, plus 3 years experience in the field of mass spectrometry to supervise new mass spectrometer laboratory.

Responsibilities will include assistance to users in solving specific problems in a wide variety of research and control applications, development and evaluation of new procedures, operation of training courses. Assist in new product developments and support marketing activities. Interest in instrumentation is a pre-requisite. We encourage publication of original work.

Excellent growth opportunity and fringe benefit program.

Forward resume to Mr. Howard M. Moran,  
Perkin-Elmer Corporation, Main Ave.,  
Norwalk, Conn. An equal opportunity employer.

**PERKIN-ELMER**

Circle No. 168 on Readers' Service Card